

# ProDetect® CA125 Rapid Test (Whole blood/Serum/Plasma)

🎯 High sensitivity, specificity and accuracy.   ✅ MDA registered.



Manufactured by:



Medical Innovation Ventures Sdn. Bhd. (MEDIVEN®)  
1st Floor, Plot 88F, Lintang Bayan Lepas 10, Bayan Lepas Industrial Park, Phase 4,  
11900 Bayan Lepas, Pulau Pinang, Malaysia.  
Email: [sales@mediven.com.co](mailto:sales@mediven.com.co) | [export@mediven.com.co](mailto:export@mediven.com.co)  
Office / Fax: +604 305 2730  
Website: [www.mediven.com.my](http://www.mediven.com.my)



## ABOUT OUR TEST KIT

A rapid and convenient immunochromatographic assay for the qualitative detection of CA125 in human whole blood, serum or plasma. It is intended for professional use as an aid in the diagnosis of ovarian cancer. This test serves only as a preliminary result. Clinical expertise and professional judgment should be sought to further evaluate the result of the test. For *in vitro* and professional use only.



- Chromatographic immunoassay.
- Detects CA125 in whole blood, serum or plasma.
- Simple test procedure.
- Rapid results in 10 minutes.



**FAST**  
Test result in 10 minutes.



**USER-FRIENDLY**  
Simple test procedure.



**RELIABLE**  
High sensitivity, specificity, and accuracy.



**FLEXIBLE**  
Able to test different types of patient specimens.

## PERFORMANCE CHARACTERISTICS

Diagnostic sensitivity and specificity

ProDetect® CA125 Rapid Test (Whole Blood/Serum/Plasma)	
Sensitivity (%)	96.00
Specificity (%)	98.00
Accuracy (%)	97.67

## ORDERING INFORMATION

Catalog No.	Product Name	No. of Tests/Kit
PR-CA125	ProDetect® CA125 Rapid Test (Whole Blood/Serum/Plasma)	25