

ProDetect® Chikungunya IgG/IgM Rapid Test

🎯 High sensitivity, specificity and accuracy. ✅ MDA registered.



Manufactured by:



Medical Innovation Ventures Sdn. Bhd. (MEDIVEN®)
1st Floor, Plot 88F, Lintang Bayan Lepas 10, Bayan Lepas Industrial Park, Phase 4,
11900 Bayan Lepas, Pulau Pinang, Malaysia.
Email: sales@mediven.com.co | export@mediven.com.co
Office / Fax: +604 305 2730
Website: www.mediven.com.my

ABOUT OUR TEST KIT

A rapid immunochromatographic test for the visual screening of Chikungunya’s IgG and IgM antibodies in human whole blood, serum or plasma as an aid in the diagnosis of Chikungunya infection. For *in vitro* and professional use only.

- Chromatographic immunoassay.
- Screening of Chikungunya’s IgG and IgM antibodies in human whole blood, serum, or plasma.
- Simple test procedure.
- Rapid results in 15 minutes.



FAST
Test result in 15 minutes.



USER-FRIENDLY
Convenient and easy test procedure.



RELIABLE
High sensitivity, specificity, and accuracy.



FLEXIBLE
Able to test different types of patient specimens.

PERFORMANCE CHARACTERISTICS

Diagnostic sensitivity and specificity

Specimen types	Antibodies	Sensitivity (%)	Specificity (%)	Accuracy (%)
Whole blood	IgM	85.0	>99.0	90.0
	IgG	88.9	>99.9	94.7
Serum	IgM	92.9	>99.9	95.4
	IgG	93.3	>99.9	96.4
Plasma	IgM	91.7	>99.9	94.1
	IgG	>99.9	>99.9	>99.9

ORDERING INFORMATION

Catalog No.	Product Name	No. of Tests/Kit
PR-CHK	ProDetect® Chikungunya IgG/IgM Rapid Test	25