

# ProDetect® Chlamydia Rapid Test (Swab/Urine)

High sensitivity, specificity and accuracy. MDA registered.



Manufactured by:



Medical Innovation Ventures Sdn. Bhd. (MEDIVEN®)  
1st Floor, Plot 88F, Lintang Bayan Lepas 10, Bayan Lepas Industrial Park, Phase 4,  
11900 Bayan Lepas, Pulau Pinang, Malaysia.  
Email: sales@mediven.com.co | export@mediven.com.co  
Office / Fax: +604 305 2730  
Website: www.mediven.com.my



## ABOUT OUR TEST KIT

A rapid immunochromatographic assay for the qualitative detection of *Chlamydia trachomatis* in female cervical swab, male urethral swab and male urine. It is intended for professional use as an aid in the diagnosis of Chlamydia infection. This test serves only a preliminary result. Clinical expertise and professional judgment should be sought to further evaluate the result of the test. For *in vitro* and professional use only.

- Chromatographic immunoassay.
- Detects Chlamydia antigen in cervical or urethral swab or male urine specimens.
- Simple test procedure.
- Rapid results in 10 minutes.



**FAST**  
Test result in 10 minutes.



**USER-FRIENDLY**  
Convenient and easy test procedure.



**RELIABLE**  
High sensitivity, specificity, and accuracy.



**FLEXIBLE**  
Able to test different types of patient specimens.

## PERFORMANCE CHARACTERISTICS

Diagnostic sensitivity and specificity

ProDetect® Chlamydia Rapid Test (Swab/Urine)	
Sensitivity (%)	92.59
Specificity (%)	97.78
Accuracy (%)	96.22

## ORDERING INFORMATION

Catalog No.	Product Name	No. of Tests/Kit
PR-CH	ProDetect® Chlamydia Rapid Test (Swab/Urine)	25