

ProDetect® Dengue Duo NS1 Ag & IgG/IgM Rapid Test

🎯 High sensitivity, specificity and accuracy. ✅ MDA registered.



Manufactured by:



Medical Innovation Ventures Sdn. Bhd. (MEDIVEN®)
1st Floor, Plot 88F, Lintang Bayan Lepas 10, Bayan Lepas Industrial Park, Phase 4,
11900 Bayan Lepas, Pulau Pinang, Malaysia.
Email: sales@mediven.com.co | export@mediven.com.co
Office / Fax: +604 305 2730
Website: www.mediven.com.my



ABOUT OUR TEST KIT

A rapid immunochromatographic test for the visual screening of dengue NS1 antigen and IgG/IgM antibodies in the human whole blood, serum or plasma as an aid in the diagnosis of dengue infection. For *in vitro* and professional use only.

- Chromatographic immunoassay.
- Detects dengue NS1 antigen and IgG/IgM antibodies in the human whole blood, serum or plasma.
- Simple test procedure.
- Rapid results in 10 minutes.



FAST
Test result in 10 minutes.



USER-FRIENDLY
Simple test procedure.



RELIABLE
High sensitivity, specificity, and accuracy.



FLEXIBLE
Able to test different types of patient specimens.

PERFORMANCE CHARACTERISTICS

Diagnostic sensitivity and specificity

	Sensitivity (%)	Specificity (%)
NS1	95.8	96.2
IgM	83.3	>99.9
IgG	>99.9	>99.9

ORDERING INFORMATION

Catalog No.	Product Name	No. of Tests/Kit
PR-DEN-3-1	ProDetect® Dengue Duo NS1 Ag & IgG/IgM Rapid Test	10
PR-DEN-3-2		25