

ProDetect® Drug of Abuse Test (Key Cup)

⊕ High sensitivity, specificity and accuracy. ✓ MDA registered.



Manufactured by:
MEDIVEN®
MEDICAL INNOVATION VENTURES

Medical Innovation Ventures Sdn. Bhd. (MEDIVEN®)
1st Floor, Plot 88F, Lintang Bayan Lepas 10, Bayan Lepas Industrial Park, Phase 4,
11900 Bayan Lepas, Pulau Pinang, Malaysia.
Email: sales@mediven.com.co | export@mediven.com.co
Office / Fax: +604 305 2730
Website: www.mediven.com.my

PRO DETECT®

ABOUT OUR TEST PANEL

A visual one-step lateral flow chromatographic immunoassay for the qualitative detection of multiple substances in human urine. This assay provides only preliminary qualitative test result. Clinical consideration and professional judgement should be applied to any drug of abuse test result, particularly when preliminary positive results are obtained. For *in vitro* diagnostic and professional use only.



FAST

Test results in 5 minutes.



USER-FRIENDLY

Simple test procedure.



RELIABLE

High sensitivity, specificity, and accuracy.




COST EFFECTIVE

Cost effective method for drug screening.

PERFORMANCE CHARACTERISTICS

Diagnostic sensitivity and specificity



ORDERING INFORMATION

Catalog No.	Product Name	No. of Test Cups/Kit
PR-KC-6	ProDetect® Drug of Abuse Test (Key Cup)	25

Tests and cut-off concentrations of drugs/drug metabolites:

Test	Cut-off (ng/ml)	Test	Cut-off (ng/ml)	Test	Cut-off (ng/ml)
Amphetamine (AMP 1000)	1,000	Methamphetamine (mAMP 500)	500	Cotinine (COT)	200
Amphetamine (AMP 500)	500	Methamphetamine (mAMP 300)	300	Synthetic Cannabinoid (K2 50)	50
Amphetamine (AMP 300)	300	Methylenedioxymethamphetamine (MDMA)	500	Synthetic Cannabinoid (K2 20)	20
Barbiturates (BAR)	300	Opiate (OPI 300, MOP, MOR)	300	Ketamine (KET)	1,000
Benzodiazepines (BZO)	300	Opiate (OPI 100)	100	Fentanyl (FEN)	200
Buprenorphine (BUP)	10	Opiate (OPI 2000)	2,000	Tramadol (TRA 100)	100
Cocaine (COC 300)	300	Oxycodone (OXY)	100	Tramadol (TRA 50)	50
Cocaine (COC 150)	150	Phencyclidine (PCP)	25	Methylenedioxypyrovalerone (MDPV)	1,000
Marijuana (THC 50)	50	Propoxyphene (PPX)	300	Ethyl Glucuronide (ETG)	300
Marijuana (THC 20)	20	Tricyclic Antidepressants (TCA)	1,000	Alcohol (ALCO)	>0.04%
Methadone (MTD)	300	2-Ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine (EDDP)	300	-	-
Methamphetamine (mAMP 1000)	1,000	6-Acetylmorphine (6-ACM)	10	-	-