

ProDetect® Flu A & B Rapid Test

High sensitivity, specificity and accuracy. MDA registered.



Manufactured by:



Medical Innovation Ventures Sdn. Bhd. (MEDIVEN®)
1st Floor, Plot 88F, Lintang Bayan Lepas 10, Bayan Lepas Industrial Park, Phase 4,
11900 Bayan Lepas, Pulau Pinang, Malaysia.
Email: sales@mediven.com.co | export@mediven.com.co
Office / Fax: +604 305 2730
Website: www.mediven.com.my

ABOUT OUR TEST KIT

A rapid immunochromatographic test for the qualitative screening of influenza A and influenza B viruses as an aid in the diagnosis of influenza infection. For *in vitro* and professional use only.

- Chromatographic immunoassay.
- Screen for influenza A and influenza B viruses in the nasopharyngeal swab, pharyngeal swab and nasopharyngeal aspirate.
- Simple test procedure.
- Rapid results in 15 minutes.



FAST
Test result in 15 minutes.



USER-FRIENDLY
Simple test procedure.



RELIABLE
High sensitivity, specificity, and accuracy.



FLEXIBLE
Able to test different types of patient specimens.

PERFORMANCE CHARACTERISTICS

Diagnostic sensitivity and specificity

Influenza type	Specimen type	Sensitivity (%)	Specificity (%)	Accuracy (%)
Influenza A	Nasopharyngeal swab	99.0	98.9	98.9
	Pharyngeal swab	95.1	99.3	98.1
	Nasopharyngeal aspirate	>99.9	99.2	99.3
Influenza B	Nasopharyngeal swab	97.7	99.0	98.6
	Pharyngeal swab	94.2	99.4	97.8
	Nasopharyngeal aspirate	97.9	99.4	98.8

ORDERING INFORMATION

Catalog No.	Product Name	No. of Tests/Kit
PR-FLU	ProDetect® Flu A & B Rapid Test	25