

ProDetect® HBsAg Rapid Test (Cassette)

🎯 High sensitivity, specificity and accuracy. ✔️ MDA registered.



Manufactured by:



Medical Innovation Ventures Sdn. Bhd. (MEDIVEN®)
1st Floor, Plot 88F, Lintang Bayan Lepas 10, Bayan Lepas Industrial Park, Phase 4,
11900 Bayan Lepas, Pulau Pinang, Malaysia.
Email: sales@mediven.com.co | export@mediven.com.co
Office / Fax: +604 305 2730
Website: www.mediven.com.my



ABOUT OUR TEST KIT

A rapid and convenient chromatographic immunoassay for the qualitative detection of hepatitis B surface antigen (HBsAg) in whole blood, serum or plasma specimens. It is intended for professional use as an aid in the diagnosis of hepatitis B. The test serves only as a preliminary result. Clinical expertise and professional judgment should be sought to further evaluate the result of the test. For *in vitro* and professional use only.

- Chromatographic immunoassay.
- Detects HBsAg in whole blood, serum or plasma.
- Simple test procedure.
- Rapid results in 15 minutes.



FAST
Test result in 15 minutes.



USER-FRIENDLY
Simple test procedure.



RELIABLE
High sensitivity, specificity, and accuracy.



FLEXIBLE
Able to test different types of patient specimens.

PERFORMANCE CHARACTERISTICS

Diagnostic sensitivity and specificity

	Sensitivity (%)	Specificity (%)	Accuracy (%)
Whole blood	>99.00	99.72	99.74
Serum	>99.90	99.50	99.64
Plasma	>99.00	99.72	99.78

ORDERING INFORMATION

Catalog No.	Product Name	No. of Tests/Kit
PR-HBsAg	ProDetect® HBsAg Rapid Test (Cassette)	25