

ProDetect[®] Rota/Adenovirus Rapid Test

🎯 High sensitivity, specificity and accuracy. ✅ MDA registered.



Manufactured by:



Medical Innovation Ventures Sdn. Bhd. (MEDIVEN[®])
1st Floor, Plot 88F, Lintang Bayan Lepas 10, Bayan Lepas Industrial Park, Phase 4,
11900 Bayan Lepas, Pulau Pinang, Malaysia.
Email: sales@mediven.com.co | export@mediven.com.co
Office / Fax: +604 305 2730
Website: www.mediven.com.my

ABOUT OUR TEST KIT

A rapid immunochromatographic test for the qualitative screening and differentiation of rotavirus and adenovirus in human fecal samples to aid in the diagnosis of rotavirus and adenovirus infection. For *in vitro* and professional use only.



- Immunochromatographic test.
- Detects antigens of rotavirus and adenovirus in a test cassette.
- Save time and convenient.
- Simple test procedure.
- Rapid results in 10 minutes.



FAST
Test result in 10 minutes.



USER-FRIENDLY
Simple test procedure.



RELIABLE
High sensitivity, specificity, and accuracy.



NON-INVASIVE
Test with human fecal samples.

PERFORMANCE CHARACTERISTICS

Diagnostic sensitivity and specificity

	Sensitivity (%)	Specificity (%)	Accuracy (%)
Rotavirus	97.4	97.2	97.3
Adenovirus	95.4	97.8	97.0

ORDERING INFORMATION

Catalog No.	Product Name	No. of Tests/Kit
PR-GAST-3	ProDetect [®] Rota/Adenovirus Rapid Test	25