

ProDetect® *Salmonella typhi* Antigen Rapid Test (Feces)

🎯 High sensitivity, specificity and accuracy. ✅ MDA registered.



ABOUT OUR TEST KIT

A rapid immunochromatographic assay for the qualitative detection of *Salmonella typhi* antigens in human feces specimens to aid in the diagnosis of *Salmonella typhi* infection. This test serves only a preliminary result. Clinical expertise and professional judgment should be sought to further evaluate the result of the test. For *in vitro* and professional use only.



- Chromatographic immunoassay.
- Detects *Salmonella typhi* antigens in human feces specimens.
- Simple test procedure.
- Rapid results in 5 minutes.



FAST
Test result in 5 minutes.



USER-FRIENDLY
Simple test procedure.



RELIABLE
High sensitivity, specificity, and accuracy.



NON-INVASIVE
Test with human feces specimens.

PERFORMANCE CHARACTERISTICS

Diagnostic sensitivity and specificity

ProDetect® <i>Salmonella typhi</i> Antigen Rapid Test (Feces)	
Sensitivity (%)	98.08
Specificity (%)	99.19
Accuracy (%)	98.86

ORDERING INFORMATION

Catalog No.	Product Name	No. of Tests/Kit
PR-STF	ProDetect® <i>Salmonella typhi</i> Antigen Rapid Test (Feces)	25

Manufactured by:



Medical Innovation Ventures Sdn. Bhd. (MEDIVEN®)
1st Floor, Plot 88F, Lintang Bayan Lepas 10, Bayan Lepas Industrial Park, Phase 4,
11900 Bayan Lepas, Pulau Pinang, Malaysia.
Email: sales@mediven.com.co | export@mediven.com.co
Office / Fax: +604 305 2730
Website: www.mediven.com.my