

ProDetect® Syphilis Rapid Test

🎯 High sensitivity, specificity and accuracy. ✔️ MDA registered.



Manufactured by:



Medical Innovation Ventures Sdn. Bhd. (MEDIVEN®)
1st Floor, Plot 88F, Lintang Bayan Lepas 10, Bayan Lepas Industrial Park, Phase 4,
11900 Bayan Lepas, Pulau Pinang, Malaysia.
Email: sales@mediven.com.co | export@mediven.com.co
Office / Fax: +604 305 2730
Website: www.mediven.com.my

ABOUT OUR TEST KIT

A rapid and convenient immunochromatographic assay used for the qualitative detection of antibodies against *Treponema pallidum* in whole blood, serum or plasma samples. It is intended for professional use as an aid in the diagnosis of Syphilis. This assay provides only as a preliminary result. Clinical expertise and professional judgment should be sought to further evaluate the result of the test. For *in vitro* and professional use only.

- Chromatographic immunoassay.
- Detects antibodies against *Treponema pallidum* in whole blood, serum or plasma.
- Simple test procedure.
- Rapid results in 5 minutes.



FAST

Test result in 5 minutes.



USER-FRIENDLY

Simple test procedure.



RELIABLE

High sensitivity, specificity, and accuracy.



FLEXIBLE

Able to test different types of patient specimens.



PERFORMANCE CHARACTERISTICS

Diagnostic sensitivity and specificity

ProDetect® Syphilis Rapid Test

Sensitivity (%)	>99.9
Specificity (%)	99.7
Accuracy (%)	99.8

ORDERING INFORMATION

Catalog No.	Product Name	No. of Tests/Kit
PR-SYP	ProDetect® Syphilis Rapid Test	25